



Original InnoVet Parts Listing

Original InnoVet Part Description	Part Number	Dealer COST US\$
revised Aug 2007 rf		
Pricing subject to change without notice.		
Air Cylinder, Damper for hinged generator tilt	W369 (Bimba 125)	\$85.88
cable, F377 to C118 Timer PCBA	E948	\$5.68
Cable, kVp display	E950	\$45.00
Cable, mAs display	E951	\$45.00
Curtain Strip, OlderTbST W200	W215	\$31.50
Foot Loop	K090	\$12.00
Grid Cabinet	W048	
Knob with pointer (Time, mA)	W344	\$15.00
Knob without pointer (kV)	A312	\$12.00
kVp Tap Switch (8 position)	02255-000	\$172.00
LED replacement kit	00312-000	\$54.60
mA Selector Switch assy	W354	\$93.00
mA selector switch cable	Supplied with W354	
Manual - (See MASTERS = Not in OBSOLETE)	K186	\$45.00
On/Off Circuit Breaker	A365	\$170.00
pcb, kVp/mAs Interlock	F377	\$490.00
pcb, Motor Start (rotor control)	A455	\$292.50
pcb, Timer Driver	C118	\$465.00
Relay, RE-1 / RE-2	HH91	\$45.50
Relay socket assy for RE-2 (exposure)	W351	\$25.00
Relay socket for RE-1 (boost)	HH81	\$5.00
Remote exposure switch (wall mount)	W413	\$128.00
Resistor, Filament 100 Ohm, 175 W RX/RXS	A347	\$24.00
Resistor, Filament 25 Ohm, 175 W RSCC	A345	\$24.00
Rotor Cap 25uF, 330V	HAB18	\$26.00
SCR Kit (includes Instructions)	00424-000	\$64.75
SCR block ^ Please order with instructions ^	00194-000	\$68.00
Sliders, Front Grid Cabinet	00673-000	\$10.00
Snubber/Suppressor (Tap Switch)	W440	\$25.00
Stabilizer cap, C2, 660V 2.5uF	A343	\$20.00
Stator Capacitor, C1, 330V 25uF(formerly A351)	HAB18	\$26.00
Switch for RemoteAssy (Cap part no. J783)	C414	\$24.84
switch button, (cap)	J783	\$2.07
Time encoder EPROM (for F377 Board)	W341	\$40.00
Time Selector Switch assy	J976	\$88.50
Time Selector Switch Cable	comes with J976	
Transformer, Space Charge	L926	\$90.00
Transformer, Stabilizer, Filament	A342	\$220.00
TS Capacitor, C2	A343	\$10.63
Two position prep/exposure foot switch	A981	\$162.50